



dNTP Mix (10 mM each)/(2.5 mM each)

NB-54-0124

NB-54-0125



dNTP Mix (10 mM each)/(2.5 mM each)

cat# NB-54-0124

cat# NB-54-0125

Introduction

The dNTP Mix (PCR Grade) is a equimolar mixture of high purity dATP, dTTP, dCTP, and dGTP, which can be directly used for PCR, reverse transcription (RT), and one-step RT-PCR.

Contents of Kits

Components	NB-54-0124-01	NB-54-0124-02	NB-54-0125-01	NB-54-0125-01
dNTP Mix (10 mM each)	1 ml	5 ml	-----	-----
dNTP Mix (2,5 mM each)	-----	-----	1 ml	5 ml

* PH 7.0-9.0.

Storage

All components should be stored at -20°C.

Quality Control

Purity > 99% (HPLC).

No exonuclease, endonuclease, and RNase activities.

Validated for general PCR, high-fidelity PCR, RT-PCR, and Real-time PCR.

Protocol

This mix can be used directly without dilution.

1. For a 50 µl PCR reaction system: add 1 µl of NB-54-0124 or 4 µl of NB-54-0125
2. For a 20 µl RT-PCR reaction system: add 1 µl of NB-54-0124 or 4 µl of NB-54-0125