

## Proteus 1-step Batch Midi Plus Spin Column Pack

#Cat: NB-45-00058-1


Size: 8pc

### Materials Supplied in the Kit:

- Proteus spin column (20 ml capacity in a swing bucket rotor)
- Two caps:
  - Clear spin push cap for all centrifugation steps
  - **Yellow** screw cap for the batch incubation step only
- 50 ml centrifuge tubes

### Additional Materials Required:

- 0.2 µm syringe filters for clarification
- 50 ml centrifuge tubes
- A bench-top centrifuge with swing bucket rotor capable of handling 50ml centrifuge tubes. (The preferred rotor is a swing bucket rotor)
- Quartz cuvettes for UV absorbance measurements
- UV/VIS spectrophotometer



### Recommended Protocol

The following spin speeds and times are appropriate for a 0.25 – 1 ml resin bed volume. Spin times foreach of the following steps may increase with larger bed volumes.

#### Pre-equilibration



1. Pipette the appropriate resin slurry into the batch incubation chamber of the spin columnbarrel. Spin the resin at 750 x g for 5 min. **This step is critical to ensure that all ethanol is removed from the resin. Many resins are stored in 20-30% ethanol. N.B. Ethanol does interfere with sealing properties of the Self Seal™ membrane technology.**



2. Pre-equilibrate the Midi spin column with 15 ml equilibration buffer by centrifuging the spin column at 750 x g for 5 min. It is **critical** that you repeat this step one more time with a

further 15 ml fresh equilibration buffer to remove any residual ethanol.

**NOTE:** If using one spin column, ensure that the spin column is counterbalanced with a unit of equal weight (adjusted with distilled water).

### Clarification of sample

3. Pre-filter the sample through a final 0.2 µm filter (e.g. syringe filter).

**NOTE:** As with all forms of chromatography, it is critical that the sample is filtered to a final 0.2 µm **immediately** before loading it on the spin column. Optimal performance of these devices will depend on these instructions being rigorously followed.

### Sample loading

4. Transfer the spin column barrel to a fresh 50 ml centrifuge tube and load the required volume of filtered sample. The maximum sample volume is 20 ml. Tightly screw the **YELLOW** batch incubation cap and invert 2-3 times to mix the sample and the resin. Place the column on a standard tube roller or rotator and mix for 1-3 hours.
5. After batch incubation, replace the yellow cap with the CLEAR spin push cap. Centrifuge the column at 750 x g for up to 10 min and collect the flow through.

**NOTE:** If using one spin column, ensure that the spin column is counterbalanced with a unit of equal weight (adjusted with distilled water).

### PURIFIED SAMPLE

6. Wash off any unbound protein with 20 ml binding buffer at 750 x g for 5 min. Repeat this step, if necessary, to ensure that all the unbound protein has been removed eg A280 < 0.1. Transfer the spin column into a fresh collection tube and then elute the target protein with up to 1-10 ml elution buffer by centrifuging the spin column at 750 x g for 5 min. The eluate contains the target protein and is now ready for further downstream analyses.