

Cat.No. SS-830  
Batch S15-11-I  
Quantity 50 µl

**Polyclonal antibody to****Somatostatin Receptor SST3 (SSTR3) Mouse / Rat**

Host: Rabbit

Immunogen: Peptide of rat and mouse SSTR3 Receptor coupled to KLH

Presentation: Whole Antiserum, liquid

Reconstitution: Reconstitute in an appropriate volume of distilled water and store in aliquots at -20°C or colder.

Specificity: This antiserum recognizes the C-terminus of rat and mouse SSTR3 receptor. Specificity of antibody binding was accomplished by blocking with its cognate peptide.

Application: Immunolocalization

Titer: Working dilution has been determined to be 1:3000 – 1:4000 in Immunohistochemistry.

Storage: Aliquot and store at -20°C or colder. Avoid repeated freeze and thawing cycles.

