

NeoBiotech

Product Description:	GHRH receptor antibody, human, C-term. Growth hormone releasing hormone receptor G Protein-coupled receptor Polyclonal antibody raised against the C-terminus of the human GHRH receptor
Product-No.:	NB-49-011-50ul
Immunogen:	residues 403-422 coupled to KLH (NP_000814) KLH- RTRAKWTTSPRSAAKVLTSM
Host:	Rabbit
Form:	serum with 0.02% NaN ₃ as preservative
Quantity:	50 µl



GHRH receptor immunoreactivity in human ovarian carcinoma.

Applications:	ICH (p), IC, W W = Western blot (1:5000) IC = Immunocytochemistry (1:5000) ICH (p) = Immunohistochemistry of formalin-fixed, paraffin-embedded tissues (1:2500)
Species cross-reactivity:	human
Sensitivity:	50 ng
Storage:	4° C short term -20° C long term

Neo-Biotech

74 rue des Suisses – 92000 Nanterre