

Product Description: CXC Chemokine Receptor 4 (CXCR4), Human, Rat, Mouse  
anti-human, rat, mouse CXCR4 antibody raised against the carboxy-terminus of the human, rat and mouse CXCR4 receptor.

Category: Cells and Surface Antigens

Order number: CXC-4000-RM

Host: Rb

Description: Rabbit IgG


Application: IHC (p), IC, WB, IP  
W = Western blot (1:1000)  
IC = Immunocytochemistry (1:500)  
IHC (p) = Immunohistochemistry of formalin-fixed, paraffin-embedded tissues (1:100)

Reactivity: Hu, Ms, Rt

Quantity: 500 µl

Immunogen: amino acid sequence 338-359  
(KGKRGHSSVSTESESSFHSS)  
of human, rat and mouse CXCR4 receptor coupled to KLH

Storage: 4 °C



CXC4 immunoreactivity in a human cervical carcinoma.