

NeoBiotech

Product Description: VPAC2 antibody, human, C-term.
VIPR2
Vasoactive intestinal peptide receptor 2
G Protein-coupled receptor
Polyclonal antibody raised against the
C-terminus of the human VPAC2 receptor


Product-No.: NB-49-013-50ul

Immunogen: residues 419-438 coupled to KLH (NP_003373)
KLH- LQFHRGSRAQSFLQTETSVI

Host: Rabbit

Form: serum with 0.02% NaN₃ as preservative

Quantity: 50 µl



VPAC2 receptor immunoreactivity in mucosa of human small intestine.

Applications: ICH (p), IC, W
W = Western blot (1:5000)
IC = Immunocytochemistry (1:5000)
ICH (p) = Immunohistochemistry of formalin-fixed, paraffin-embedded tissues (1:2500)

Species cross-reactivity: human

Sensitivity: 50 ng

Storage: 4°C short term
-20°C long term

Neo-Biotech

74 rue des Suisses – 92000 Nanterre