


## Tris(2-carboxyethyl)phosphine hydrochloride

[51805-45-9]

Cat # NB-42-02230



### Product Information

CAS : 51805-45-9

Sizes 10g, 25g

Molecular Weight 286.65 g/mol

Melting Point 177-178 °C

Chemical Hazard Corrosive

Hazard Statements H314: Causes severe skin burns and eye damage.

Precautionary Statements  
P260: Do not breathe.  
P271: Use only outdoors or in a well-ventilated area.  
P280: Wear protection equipment hazard.

*For reference only  
For Research Use Only. Not for Diagnostic or Therapeutic Use.*