



Anti-RXFP2 antibody (90-170 Internal)

Cat# NB-22-62344-1

Cat# NB-22-62344-50

Cat# NB-22-62344-500

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HLA-A, B, C is suitable for use in Flow Cytometry, Western Blot, Immunohistochemistry and Immunoprecipitation research applications
Applications	FC, WB, IHC, IP
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	monoclonal
Clone ID	W6/32
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG2ak
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product

TARGET INFORMATION

Gene ID	3106
Gene Symbol	HLA-B
Uniprot ID	HLAB_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic other medical purposes.