

Orbital Shaker N1000A/N1000B

#Cat: NB-12-0015A (N1000A)

#Cat: NB-12-0015B (N1000B)



Product Introduction

N1000A and N1000B orbital shaker are powerful variable speed shaker which provides efficient orbital motion. Speed and time are under microprocessor control. N1000A / N1000B shaker can be applicable in different labs: microbiology, chemistry, immunology, biochemistry and molecular biology.

Features

1. Timer function, set the working time within 1 min-100 h
2. Normal operation at 4~45 °C, safe and convenient
3. A variety of pallets are available for different mixing operation

Specification

Model	NB-12-0015B	NB-12-0015A
Shaking Speed	50-300 rpm	50-250 rpm
Orbit	10 mm	20 mm
Timer range	1 min~99 h59 min	1 min~99 h59 min
Load capacity	See below table	See below table
Max. capacity	3.0 kg	2.5 kg
Motor	DC brushless	DC brushless
Operating ambient	4~45 °C	4~45 °C
Dimension (mm)	280×270×110 mm	280×270×110 mm

Ordering Information

Code	Description
NB-12-0015B	Orbital Shaker N1000B @50~300rpm
NB-12-0015A	Orbital Shaker N1000A @50~250rpm

N1000A / N1000B Platforms



Suitable Blocks

NB-12-0015-01	Flat Platform for N1000A & B
NB-12-0015-02	Platform 100ml for N1000A & B
NB-12-0015-03	Platform 250ml for N1000A & B
NB-12-0015-04	Universal Platform for N1000A & B
NB-12-0015-05	Platform with springs for N1000A & B