

NeoBiotech

Product Description: CCKB antibody, human, C-term.
CCK2
Cholecystokinin type B receptor
Cholecystokinin type 1 receptor
G Protein-coupled receptor
Polyclonal antibody raised against the
C-terminus of the human CCKB receptor

Product-No.: NB-49-004-50ul

Immunogen: residues 428-447 coupled to KLH (NP_795344)
KLH- psiaslsrslsyttistlgpg

Host: Rabbit

Form: serum with 0.02% NaN as preservative


Quantity: 50 µl

Applications: IC, W
W = Western blot (1:5000)
IC = Immunocytochemistry (1:5000)

Species cross-reactivity: human

Sensitivity: 50 ng

Storage: 4 °C short term
-20 °C long term



Neo-Biotech

74 rue des Suisses – 92000 Nanterre