


## Plasmid To Express Inntag 5

Cat # NB-19-0048

Inntag #	pdb_id	Description	prot_gi	Organism	Length (aa)
5	1WHP	Allergen Phl p 2	157834206	<i>Phleum pratense</i>	96



### CLONING INFORMATION

Cloning method Restriction Enzyme:

5' cloning site: EcoRI (not destroyed)

3' cloning site: SalI (not destroyed)

5' sequencing primer: GTTGTG TGGAATTGTGAGCG

### GROWTH IN BACTERIA

Bacterial Resistance(s): Ampicillin

Growth Temperature: 37°C

Growth Strain(s): DH5alpha

Copy number: High Copy

### Reference:

Inntags: small self-structured epitopes for innocuous protein tagging. Georgieva MV, Yahya G, Codo L, Ortiz R, Teixido L, Claros J, Jara R, Jara M, Iborra A, Gelpi JL, Gallego C, Orozco M, Aldea M. Nat Methods. 2015 Aug 31. doi: 10.1038/nmeth.3556. 10.1038/nmeth.3556 PubMed 26322837

For Research Use Only. Not for Diagnostic or Therapeutic Use.