

# Product Data Sheet



74, rue des Suisses

92000 Nanterre, FRANCE Tel: (33) 9 77 40 09  
09 | Fax: (33) 9 77 40 10 11  
Email: [info@neo-biotech.com](mailto:info@neo-biotech.com)  
<https://www.neo-biotech.com>

<b>Product</b>	<b>Sheep Whole Blood (RBC), defibrinated</b>
<b>Product No</b>	NB-43-0795-10
<b>Lot</b>	31-112-060KW02/21
<b>Date</b>	12/01/2021
<b>Volume</b>	1 x 10 mL
<b>Presentation</b>	1 lot of 10 mL of sheep whole blood, defibrinated, fresh, pooled from 1 female and 1 male animal
<b>Species of Origin</b>	Sheep
<b>Strain</b>	Merino Land Race
<b>Sex</b>	female, male
<b>Remarks</b>	non medicated, non immunized, non pregnant female sheep
<b>Anticoagulant</b>	non
<b>Sterility</b>	aseptic
<b>Buffer</b>	non
<b>Preservative</b>	non
<b>Storage Conditions</b>	+2 - +8 °C
<b>Expiry Date</b>	16/01/2021
<b>Disclaimer</b>	<p>THIS MATERIAL SHOULD BE HANDLED AS IF CAPABLE OF TRANSMITTING INFECTIOUS AGENTS. PLEASE USE UNIVERSAL PRECAUTIONS.</p> <p>THIS PRODUCT IS BEING SOLD FOR RESEARCH AND, OR MANUFACTURING PURPOSES ONLY. IT IS NOT TO BE USED IN HUMANS OR ANIMALS AND FURTHER MANUFACTURING THAT WILL RESULT IN A FINISHED PRODUCT THAT CONTAINS VIABLE LEUKOCYTES IS PROHIBITED. FOR IN VITRO USE ONLY. THE USER ASSUMES ALL RESPONSIBILITY FOR ITS USAGE AND DISPOSAL, IN ACCORDANCE WITH ALL REGULATIONS.</p>
<b>Country of Origin</b>	<b>Germany</b>