

## Mini Centrifuge N600 @4,000 & 6,000rpm (Including 2 quick release interchangeable rotors for 8x1.5/2.0ml tubes & 4 place PCR Strip)

#Cat: NB-12-0045

N600 is an improved model of N500A. It has large capacity and can place 8-well 1.5 mL/2.0 mL tubes or 4 rows of 8x0.2 mL PCR strip tubes. Two speeds are interchangeable, suitable for the two rotors.

### Features

- ▶ Compact and streamlined design, small size
- ▶ Mute technology, ultra-low noise operation
- ▶ Using micro DC motor and carbon fiber rotor
- ▶ 8-well rotor and 4x8x0.2 mL large-capacity rotor design
- ▶ Quick release rotor tub, makes rotor exchange easier
- ▶ Two speed options: 4,000 & 6,000 rpm



Two Speed Options

### N600 Mini Centrifuge

#### Rotor

- ▶ 8x1.5 mL/2.0 mL tube rotor
- ▶ 4x8x0.2 mL PCR strip rotor



### Specification

Model	Mini-6KC
Max. RCF	2,000xg (6,000 rpm) 1,200xg (4,000 rpm)
Capacity	8x1.5 mL/2.0 mL tubes 4x8x0.2 mL PCR strip tubes 8 or 32x0.2 mL tubes 8x0.5 mL tubes
Speed	6,000 rpm & 4,000 rpm
Noise	<50 dB

### Ordering Information

Code	Description
NB-12-0045	N600 Mini centrifuge
NB-12-0045-04	8x1.5 mL/2.0 mL rotor
NB-12-0045-03	4x8x0.2 mL strip rotor
NB-12-0045-01	Adapters for 0.5 mL tubes
NB-12-0045-02	Adapters for 0.2 mL tubes