

Neo Biotech

Product Description: NMBR antibody, human, C-term.
Neuromedin B receptor
Neuromedin B-preferring bombesin receptor
Bombesin receptor
G Protein-coupled Receptor
Polyclonal antibody raised against the carboxy-terminus of the human NMBR receptor.

Product-No.: NB-49-009-50ul

Immunogen: residues 342-361 coupled to KHL (CAH70473)
KLH-nmvtnsvllnghsmkqema1

Host: Rabbit

Form: Antiserum

Quantity: 50 µl

Applications: IC, W
W = Western blot (1 µg/ml)
IC = Immunocytochemistry (1 µg/ml)

Species cross-reactivity: human

Storage: 4 °C short term
-20 °C long term