

Cat.No. SS-800  
Batch S15-11-I  
Quantity 50 µl

**Polyclonal antibody to  
Somatostatin Receptor SST2A (SSTR2A)**

Host: Rabbit  
Immunogen: Peptide coupled to KLH  
Presentation: Whole Antiserum, liquid  
Specificity: This antiserum recognizes the C-terminus of human, rat and mouse SSTR2A receptor. Specificity of antibody binding was accomplished by blocking with its cognate peptide.  
Application: Immunohistochemistry  
Titer: Working dilution has been determined to be 1:3000 – 1:4000 in Immunohistochemistry.  
Storage: Aliquot and store at -20°C or colder. Avoid repeated freezing and thawing cycles.

