



DATE: 01SEP2023

PRODUCT DATA SHEET

CYNOMOLGUS MONKEY (EU) SERUM (POOLED GENDER)

Product Number: NB - 53 - 0023 - 50

Lot Number(s): CYN351414
Collection Date: 03JUN2022
Expiration Date: 03JUN2027
Storage Conditions: -20°C OR COLDER
Anticoagulant: NONE

PROJECT DETAILS

Purchase Order:	153895
Sales Order Number:	S82023001454
Total Volume Shipped:	50 ML
Presentation:	1 LOTS OF 50ML (1 X 50ML) OF CYNOMOLGUS MONKEY (EU) SERUM, POOLED FROM 2 MALE & 2 FEMALE ANIMALS
Remarks:	NO FREEZE/THAW CYCLES Samples originate from Mauritian Cynomolgus macaques that are naturally B.virus free. All animals going into colonies are tested free of Herpes B. virus, SRV, SIV and STLV. All colonies are measles vaccinated. Mauritius is free of malaria and rabies.
Strain:	MACACA FASCICULARIS
Filtration:	None

LOT #	UID'S
CYN351414	M:SB211, YS16D F:AR667, LP81A

BIOHAZARD INFO: THIS MATERIAL SHOULD BE HANDLED AS IF CAPABLE OF TRANSMITTING INFECTIOUS AGENTS. PLEASE USE UNIVERSAL PRECAUTIONS. THESE MONKEYS HAVE BEEN INDIVIDUALLY TESTED FOR B-VIRUS (HERPES VIRUS SIMIAE) AND FOUND TO BE NEGATIVE. THIS LOT HAS BEEN PRODUCED FROM THOSE MONKEYS ONLY.

THIS PRODUCT IS BEING SOLD FOR RESEARCH AND/OR MANUFACTURING PURPOSES ONLY. THE BIOLOGICAL SAMPLES SUPPLIED BY NEO BIOTECH, OR ANY MATERIAL ISOLATED FROM THE SAMPLES, ARE FOR IN-VITRO RESEARCH USE ONLY AND ARE NOT TO BE USED AS A SOURCE OF MATERIAL FOR CLINICAL THERAPIES. HUMAN MATERIAL MAY BE USED IN VIVO IN ANIMALS. THE USER ASSUMES ALL RESPONSIBILITY FOR ITS USAGE AND DISPOSAL, IN ACCORDANCE WITH ALL REGULATIONS.

NEO BIOTECH WARRANTS THAT AT THE TIME OF SHIPMENT THE PRODUCT IS FREE FROM DEFECTS IN MATERIAL AND WORKMANSHIP AND CONFORMS TO SPECIFICATIONS WHICH ACCOMPANY THE PRODUCT. NEO BIOTECH MAKES NO OTHER WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE OR PERFORMANCE.