

Mouse Anti -Human IgG2 antibody {FITC}

#Cat: NB-22-31353-100 Size: 100µl

General Information

Product Type	Secondary antibodies
Short Description	FITC-conjugated mouse monoclonal anti-Human IgG2 secondary antibody. For use in most research applications.
Applications	WB, IHC-P, IF, ICC, ELISA
Host/Source	Mouse
Reactivity	Human

Product Properties

Clonality	Monoclonal
Clone ID	
Concentration	4-5 mg/ml
Conjugation	FITC
Purification	The antibody was isolated from ascitic fluid by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	WB 1:5000-20000 ELISA 1:10000-20000 IHC 1:5000-10000 IF 1:5000-10000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide
Isotype	IgG2
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

Target Information

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.