

Normal Rabbit serum blocking buffer

#Cat: NB-23-00189-1

Size: 100ml

#Cat: NB-23-00189-2

Size: 18ml

Intended Use:

Normal Rabbit Serum Blocking Buffer is used in immunohistochemical staining procedures to block non-specific staining when the host species of the secondary antibody is Rabbit.

Formulation:

Solution contains 0.05% NaN₃.

Recommended Protocol

1. Apply 100uL or enough volume of Normal Rabbit Serum Blocking Buffer on tissue before adding primary antibody and incubate for 10 minutes at room temperature.
2. Blot off Normal Rabbit Serum Blocking Buffer Solution (Do not rinse) and add primary antibody.

Storage

Store at 2-8°C.

Precaution:

Wear gloves and take other necessary laboratory safety procedure

For Research Use Only. Not for Diagnostic or Therapeutic Use.