

Mini Vortex

#Cat: NB-12-0039 Size: 1pcs
 #Cat: NB-12-0039-01 Size: 1pcs
 #Cat: NB-12-0039-02 Size: 1pcs

Product Introduction :

The Mini vortex is the ideal vortex mixer for every lab with its compact size, ease of use, and its economic price point. The 2,500 rpm fixed speed ensures there is no sacrificing any mixing power. It requires little working space and is easy-to-use.

Features :

1. Touch activated vortex mixer
2. Stable compact design
3. 2,500 rpm fixed vortexing speed
4. Powerful vortexing for tubes up to 15 ml
5. A power adapter and a battery case are included



Model	Mini Vortex
Cat #	NB-12-0039
Working style	Continous / Short-Spin
Shaking movement	Orbital, 3 mm
Speed Range	800 - 2500 rpm
Dimensions (mm)	74 x 91 x 84



Accessories, included in the Mini Vortex	Cat #
DC6V Power Adapter	NB-12-0039-01
Battery Box for 4 pcs AA batteries	NB-12-0039-02

