

## Sample concentration N600

#Cat: NB-12-0024A

#Cat: NB-12-0024B

#Cat: NB-12-0024C

The instrument consists of a base and stand assembly, sample holder and gas distribution system.

### Features:

1. Blowing a large number of samples at one time
2. Easy operation, safe and reliable to use, amount of airflow can be adjusted
3. High precision and wide range of temperature control



N600A



N600B



N600C

### Specifications:

Model	N600A	N600B	N600C
Cat#	NB-12-0024A	NB-12-0024B	NB-12-0024C
Temperature control range	RT.+5°C~150°C		
Temp. control accuracy @40~100°C	±0.5°C		
Temp. control accuracy @100-150°C	±1°C		
Temp. accuracy for blocks @40°C	±0.3°C		
Display Accuracy	0.1°C		
Heating Time (from 40-150°C)	≤30min		
Standard Blocks	1	1	2
Needle Length	150 mm	80 mm	150 mm
Nitrogen Consumption	0-15L/min		
Dimension (mm)	280x240x500		

**Accessories:**

NB-12-0024-01	Needles 150mm for N600A and N600C
NB-12-0024B-05	Needles 80mm for N600B
NB-12-0024B-01	Flow Control Valve for N600B
NB-12-0024B-04	Block Microplate for N600B
NB-12-0024B-03	Block 96 x 0.2 ml for N600B
NB-12-0024-02	Flow Control Valve for N600A and N600C

**Blocks:**

Code	Description	Dimension (mm)
NB-12-0009-01	49x6mm	96 x 96 x 49
NB-12-0009-02	49x7mm	96 x 96 x 49
NB-12-0009-03	25x10mm	96 x 96 x 49
NB-12-0009-04	25x12mm	96 x 96 x 49
NB-12-0009-05	25x13mm	96 x 96 x 49
NB-12-0009-06	12x15mm	96 x 96 x 49
NB-12-0009-07	16x15mm	96 x 96 x 49
NB-12-0009-08	12x16mm	96 x 96 x 49
NB-12-0009-09	16x16mm	96 x 96 x 49
NB-12-0009-10	12x19mm	96 x 96 x 49
NB-12-0009-11	16x19mm	96 x 96 x 49
NB-12-0009-12	9x20mm	96 x 96 x 49
NB-12-0009-13	9x26mm	96 x 96 x 49
NB-12-0009-14	4x28mm	96 x 96 x 49
NB-12-0009-15	4x40mm	96 x 96 x 49
NB-12-0009-16	49x0.5ml tubes	96 x 96 x 49
NB-12-0009-17	25x1.5ml tubes	96 x 96 x 49
NB-12-0009-18	25x2.0ml tubes	96 x 96 x 49
NB-12-0027-01	12 x 10mm	192 x 96 x 49
NB-12-0027-02	12 x 12mm	192 x 96 x 49
NB-12-0027-03	12 x 13mm	192 x 96 x 49
NB-12-0027-04	12 x 15mm	192 x 96 x 49
NB-12-0027-05	12 x 16mm	192 x 96 x 49
NB-12-0027-06	12 x 19mm	192 x 96 x 49
NB-12-0027-07	12 x 20mm	192 x 96 x 49

REMARK: N600B can only hold NB-12-0024B-03 and NB-12-0024B-04



Block



Valve

