

Cat.No. SS-870
Batch G VIII
Quantity 50 µl

**Polyclonal antibody to
Somatostatin Receptor, type SSTR2A**

Host: Guinea pig

Immunogen: Amino acid sequence 355-369
(ETQRTLLNGDLQTSI) of human, rat and mouse SSTR2A
Receptor; (Vanetti, M. et al. FEBS Lett.331 (1993), 260-
266) coupled to KLH

Presentation: Whole Antiserum, liquid

Specificity: This antiserum recognizes the C-terminus of human, rat
and mouse SSTR2A receptor. Specificity of antibody
binding was accomplished by blocking with its cognate
peptide.

Application: Immunolocalization of SSTR2A receptor

Titer: Working dilution has been determined to be >1/4000 in
immunocytochemistry.

